

- Begin application in a limited area. A slight tingling is normal. If the patient experiences intense burning, this means that the skin is irritable. In this case, rinse the treated area with water, stop the session and apply WiQo MOISTURIZING AND NOURISHING FACE CREAM for dry skin.
- Always massage the product on well lubricated skin to avoid excoriations (if the glove or the spatula do not glide anymore, apply more product).
- Avoid the area under the eyes (where bags form), in order to avoid unattractive swelling in this area in the days following the treatment.
- Carry out the treatment after a suitable length of time has passed since shaving or waxing.
- In order to avoid allergic reactions, use nitrile rather than latex gloves.
- Do not treat areas affected by seborrhoeic dermatitis (e.g. sides of nose).
- Do not treat the corners of the mouth if they are affected by angular cheilitis.
- Be careful when treating the inframammary fold in order to avoid the development of fissures in the event of damaged skin.
- Warning: when treating seborrhoeic or acneic skin, especially the delicate areas of the face (e.g. sides of the mouth), perform a maximum of two successive applications of PRX-T₃₃, massaging gently. Overly aggressive massaging may abrade the damaged skin and possibly cause scarring.
- If retinoic acid (isotretinoin) is used topically, suspend its use 15 days before performing the treatment with PRX-T₃₃. After completing all the PRX-T₃₃ sessions, use of retinoic acid can be resumed.



Direct line to contact dermatologists and research biologists
for scientific questions or enquiries regarding outpatient treatment:
scientificdiv.gpq@wiqo.eu

GPQ srl

Operational Headquarters: Via Flavia, 23/1 - 34148 Trieste - Italy
Tel.: +39 040 360408 - Fax: +39 040 761393 - gpq@wiqo.eu
www.wiqo.com

LOGO

DISTRIBUTOR
contacts
contacts

PRX-T₃₃®

wiQO

PROTOCOL FOR USE



REVITALISATION OF THE FACE AND NECK

Indications:

Photo-ageing, dermal atrophy in middle-aged or older patients.

Frequency of sessions:

Once a week for 1-5 sessions.

**METHOD OF USE**

- 1) Cleanse the skin with WiQomed **P SOLUTION**.
- 2) Apply a small amount of **PRX-T₃₃** on a limited part of the area to be treated and massage it with the rubber-tipped spatula or manually with the specific technique taught, being careful to keep the area lubricated with the product.
- 3) Repeat the application until the desired effect has been achieved (visible and palpable firmness). It generally takes about three applications: 1-2 applications are necessary for the treatment of the neck or for very thin and/or delicate skin, and up to 5 are needed for very thick skin.
- 4) Proceed section by section until the entire surface has been treated (avoiding the eyelids).
- 5) Remove excess gel with wet cotton wool. Have the patient rinse the skin with cold running water, thus allowing him/her to immediately feel the smoothness and firmness of the skin.



- 6) Liberally apply WiQo **MOISTURIZING AND NOURISHING FACE CREAM** for dry skin to the area treated in order to restore the hydrolipidic film, prevent dryness and soothe any temporary redness.



- 7) Prescribe post-treatment at home with WiQo **MOISTURIZING AND NOURISHING FACE CREAM** for dry skin, to be applied 2-3 times a day for the first few days and once/twice a day thereafter. Also prescribe home treatment with WiQo **SMOOTHING FACE FLUID**, to be applied once a day starting the day after **PRX-T₃₃** treatment. Remind the patient that no other cream should be used within one hour before or after the application of WiQo **SMOOTHING FACE FLUID**.

REVITALIZATION OF DÉCOLLETAGE SKIN AND OTHER AREAS OF THE BODY

Indications:

Photo-ageing, dermal atrophy in middle-aged or older patients.

Frequency of sessions:

4-6 weekly sessions.

METHOD OF USE

- 1) Apply **PRX-T₃₃** by manually massaging it into the skin with the technique taught or with the special rubber-tipped brush, down to the upper quadrants of breast.
- 2) Before the product is fully absorbed and the glove or spatula do not glide anymore, apply more product. The desired effect (visible, palpable firmness) is achieved after 3-5 applications.
- 3) Rinse with water and apply WiQo **MOISTURIZING AND NOURISHING FACE CREAM** for dry skin.
- 4) Prescribe WiQo **ELASTICIZING ANTI-DRYNESS BODY CREAM**, to be applied at home once a day every day.



STRETCH MARKS

Indications:

Recent stretch marks (pinkish or white), old stretch marks (including those associated with fractional laser treatment).

Frequency of sessions:

Once a week for 5-15 sessions.

**METHOD OF USE**

- 1) Rub **PRX-T₃₃** lengthwise into the stretch marks using the special rubber – tipped brush.
- 2) Repeat until the bottom of the stretch marks rises. This occurs on average after 3-5 applications for pinkish or light marks and 5-10 for others.
- 3) Remove excess product with water and apply WiQo **MOISTURIZING AND NOURISHING FACE CREAM** for dry skin.



RECENT DEPRESSED SCARS

Indications:

Recent depressed scars (acne, impetigo, chicken pox, etc.) which still have a pinkish base. Treatment should be started after the complete re-epithelialization of the lesion.

Frequency of sessions:

3 times a week for 4-8 weeks until the base of the scar remains raised at a check after 15 days.

METHOD OF USE

- 1) Rub the depressed base of the scar with the special rubber – tipped brush or with a cotton bud soaked in **PRX-T₃₃**.
- 2) Repeat the application until the base of the scar lifts up.
- 3) Remove excess gel with damp cotton wool and apply WiQo **MOISTURIZING AND NOURISHING FACE CREAM** for dry skin.

**General contraindications in the application of PRX-T₃₃**

- Irritated skin
- Irritable skin
- Seborrhoeic dermatitis
- Allergy to kojic acid
- Allergy to latex, if latex gloves are used
- Pregnancy and breastfeeding
- Cutaneous and/or herpetic eruptions

Possible side effects and problems:

- **Severe peeling:** this can occur after the first session of treatment on skins with a very thick corneal layer. It may last 2-3 days. Apply WiQo **MOISTURIZING AND NOURISHING FACE CREAM** for dry skin several times a day.
- **Delicate skin that becomes red at the first application:** stop after one application and recommend WiQo **MOISTURIZING AND NOURISHING FACE CREAM** for dry skin.
- **Intense burning upon application:** this can occur with patients with irritable skin. The diagnosis of "irritable skin" is symptomatological and is given when the patient complains of burning after the application of an acidic product which would be tolerated by most people. If this happens, rinse the skin immediately with cold running water and apply WiQo **MOISTURIZING AND NOURISHING FACE CREAM** for dry skin.
- **Reddening upon application:** redness will fade in a few hours with the use of WiQo **MOISTURIZING AND NOURISHING FACE CREAM** for dry skin.
- **Delayed widespread or patchy redness:** this can occur the day after treatment and is treated with cortisone cream applied for several days (in which case sun exposure should be avoided).
- **Dryness:** apply WiQo **MOISTURIZING AND NOURISHING FACE CREAM** for dry skin more frequently.
- **Allergy to kojic acid or latex gloves:** prescribe cortisone cream.
- **Contact with the conjunctiva:** wash with sterile saline solution.
- **Itchiness, with or without redness, in the days following treatment:** very rare, presumably due to an allergy to kojic acid. Prescribe a steroidal cream.
- **Delayed limited inflammation (appearing the following day or later):** redness and/or swelling of a limited area of the treated skin, presumably due to an irritation of skin that was not completely undamaged at the time of application (subclinical dermatitis, recently waxed skin, etc.). Prescribe a steroidal cream.
- **Appearance of rounded darker patches in the two days following treatment:** these are due to drying of thicker areas of skin, and will disappear through natural exfoliation in a few days.

In the event of side effects that cause persistent inflammation (erythema lasting more than one day), there is a risk of hyperpigmentation in predisposed subjects. This is not due directly to the trichloroacetic acid, but rather to the erythema. Monitor such patients to ensure that they avoid sun exposure and use a high protection sunscreen. Prescribe tyrosinase inhibitors if necessary.