

3 STEPS TO REVITALIZE THE SKIN

- painless
- needle-free
- no downtime
- with immediate results

1. PREPARE

P solution

Cleanse the skin with WiQo P solution and pat dry.



2. STIMULATE

PRX-T₃₃

Massage PRX-T₃₃ onto the skin using a specific technique until fully absorbed.



3. HYDRATE

WiQo Nourishing and Moisturizing Cream for Dry Skin

Apply WiQo Nourishing and Moisturizing Cream for Dry Skin to restore the hydrolipidic film, prevent dryness, and address any temporary redness.



Concept and art direction by Hangar Design Group

PBRPT114_V1 09/2024

Copyright (C) 2024 - WiQo S.p.A.

DIRECT LINE

with the scientific team for inquiries and complaints: scientificdivision@wiqo.eu

WiQo S.p.A.

Via Martinelli, 16
34015, Muggia (Trieste)
Italy
Ph: +39 040 360408
wiqo.com

PRX-T₃₃[®]

WIQO

NEEDLE-FREE
BIOREVITALIZATION

PROTOCOLS OF USE

- / Cutaneous elastosis of the face, neck, décolleté, and body
- / Hyperpigmentation and melasma
- / Scar atrophy (scars and stretch marks)

PIONEERS IN
NON-INVASIVE
AESTHETIC
MEDICINE



We are innovators in advanced aesthetic dermatology, formulating solutions for all skin types.

Our goal is to boost confidence in one's skin by delivering exceptional results through safe formulas and treatments for all patients.



Informational material intended exclusively for medical professionals

WiQO[®]
FORMULATORS OF SKIN CONFIDENCE

WiQO[®]
FORMULATORS OF SKIN CONFIDENCE

THE 3 PRX-T₃₃ PROTOCOLS

CUTANEOUS ELASTOSIS OF THE FACE, NECK, DÉCOLLETÉ, AND BODY:

Photoaging, atrophic dermis, skin laxity

NUMBER OF SESSIONS*: Once a week, 3 to 5 sessions for the face, and up to 4-6 sessions for body treatments.

HOME-CARE:

WiQo Nourishing and Moisturizing Cream for Dry Skin, to be applied twice a day for the first few days after treatment to prevent peeling and dryness.

WiQo Smoothing Fluid, to be applied 24 hours after treatment, once a day.

WiQo Firming Body Cream, to be applied 2-3 times a day after treatment and 1-2 times a day in the following days.

HYPERPIGMENTATION:

Post-inflammatory and melasma

NUMBER OF SESSIONS*: One session every 15 days, up to 6 sessions.

HOME-CARE:

WiQo Nourishing and Moisturizing Cream for Dry Skin, to be applied twice a day for the first few days after treatment to prevent peeling and dryness, and then once a day.

WiQo Lightening Serum, to be applied twice a day, morning and evening.

SCAR ATROPHY:

Depressed or hollow scars and stretch marks


NUMBER OF SESSIONS*: Once a week for 5-10 sessions.

HOME-CARE:

WiQo Firming Body Cream, to be applied 2-3 times a day for the first few days after treatment and 1-2 times a day in the following days.

*The number of sessions and recommended applications are guidelines at the discretion of the physician after patient evaluation.

CORRECT METHOD OF PRX-T₃₃ APPLICATION



1. To avoid allergic reactions, use nitrile gloves, not latex. Cover the tip of the WiQo spatula with a **nitrile glove**.
2. Begin application on a small, limited area.
3. Always massage the product on well-lubricated skin to avoid abrasions (if the glove or spatula no longer glides, apply more product).
4. A slight tingling sensation is normal. If the patient experiences intense burning, it indicates the skin is irritable. In such cases, rinse the area with water, stop the session, and apply **WiQo Nourishing and Moisturizing Cream for Dry Skin**.
5. When applying to acne-prone skin, especially on delicate facial areas (e.g., sides of the mouth), perform a maximum of two applications of PRX-T₃₃, massaging gently. An excessively aggressive massage can abrade damaged skin and cause scarring.
6. Avoid the under-eye area (where bags form) to prevent unwanted swelling in this area in the days following application or apply only one layer with light pressure.
7. Exercise caution when applying to the intermammary sulcus.
8. Advise the patient to use adequate sun protection after treatment.

REVITALIZE YOUR PATIENTS' SKIN IN 5 SIMPLE STEPS

STEP 1

Cleanse with **P solution** and dry thoroughly.

STEP 2

Extract PRX-T₃₃ by piercing the vial's cap with the tip of the cannula or an 18G needle. Insert the syringe into the cannula and draw out the product. Apply a few drops to the skin using the syringe without the needle.

STEP 3

Massage PRX-T₃₃ onto the skin using a specific technique until fully absorbed. Always massage on well-lubricated skin to avoid abrasions. If the glove no longer glides on the skin, apply more product.

STEP 4

Repeat the application until a visible and palpable turgor is achieved. Typically, 2-3 layers are required: 1-2 for neck treatment or sensitive/delicate skin, and up to 5 for very thick skin.

In case of scar atrophy, the correct application technique is as follows: rub PRX-T₃₃ along the scar (or stretch mark) lengthwise using the WiQo spatula.

STEP 5

Remove any residual product with a water-soaked cloth. Apply **WiQo Nourishing and Moisturizing Cream for Dry Skin** to restore the hydrolipidic film, prevent dryness, and soothe any mild redness.

WARNINGS*

- For external use only, topical application.
- Avoid contact with eyes; in case of contact, rinse with saline solution.
- If the patient is using topical retinoic acid (isotretinoin), discontinue its use 15 days before applying PRX-T₃₃.
- Do not apply to the corners of the mouth if affected by angular cheilitis.
- Avoid application on nevi of any histological type and/or on skin with burns of various degrees, sunburn, etc.
- Treat at a reasonable distance from the last shave or waxing.
- Store PRX-T₃₃ at 2°C/8°C, even after opening the vial.

*It is recommended to follow the instructions as per the package leaflet.

GENERAL CONTRAINDICATIONS

- / Irritated and/or irritable skin
- / Seborrheic or atopic dermatitis
- / Macerated skin, open wounds, rashes, and/or herpetic eruptions
- / Epidermal and dermal tumor pathologies or other neoplastic conditions with possible local involvement (e.g., cutaneous lymphoma)
- / Allergy to the product components
- / Pregnancy and breastfeeding

POSSIBLE SIDE EFFECTS

- / Exfoliation (which may be intense)
- / Redness (immediate or delayed)
- / Burning sensation during application (which may be intense and immediate)
- / Dryness
- / Anaphylaxis
- / Itching
- / Persistent inflammation
- / Appearance of dark, round patches within two days after application: these are areas of thicker skin that will dry out and disappear through natural exfoliation within a few days.

USEFUL INFORMATION

In case of side effects that cause persistent inflammation (erythema lasting more than one day), it is important to be aware of the risk of hyperpigmentation in predisposed individuals. This is not directly caused by trichloroacetic acid, but rather by the erythema. We recommend monitoring such patients to ensure they avoid sun exposure and use high-protection sunscreen. If necessary, treat the inflammation accordingly.