

Warnings:

- Begin the application in a limited area. A slight tingling is normal. If the patient experiences intense burning, it means that the skin is irritable. In this case, rinse the treated area with water, stop the session and apply WiQo Restoring Cream for External Delicate Areas.
- Always apply the product on well lubricated skin **to avoid excoriations**.
- Carry out the treatment after a suitable length of time has passed since shaving or waxing.
- In order to avoid allergic reactions, use nitrile rather than latex gloves.
- Do not treat areas affected by seborrhoeic dermatitis.
- Do not treat the inguinal fold.
- Warning: in case of contact with the eyes, rinse thoroughly with saline solution.

Possible side effects and problems:

- **Itchiness, with or without redness, in the days following the treatment:** very rare, presumably due to an allergy to kojic acid. Prescribe a steroidal cream.

WiQO^{med}
Scientific Division

A direct line to contact dermatologists and research biologists **WiQO^{med}** for scientific questions or enquiries regarding outpatient treatment

www.wiqomed.com
scientificdiv.gpq@wiqo.com

M11US0MM2_EN_V1_06/2018
Copyright © 2018 - GPQ srl
Informational material reserved for medical professionals. Reproduction prohibited.

PRX-T[®] lady

PROTOCOL FOR USE OF THE MEDICAL DEVICE



INFORMATIONAL MATERIAL RESERVED FOR MEDICAL PROFESSIONALS

WiQO^{med}



GPQ Srl - Milan
Production Laboratory and R&D
Ph.: +39 040 360408 - Fax: +39 040 761393
contact@wiqomed.com
Via Flavia, 23/1
TRIESTE - ITALY

Reproduction prohibited.

REVITALISATION AND LIGHTENING OF EXTERNAL DELICATE AREAS

Frequency of sessions:

1 session every 15 days for 5 or more sessions.

METHOD OF USE

- 1) Clean the area to be treated with WiQomed P Solution.
- 2) Perforate the cap of the vial with a needle (18g).
- 3) Insert the syringe on the needle and extract the desired quantity.
- 4) Protect the mucous membranes with WiQo Restoring Cream for External Delicate Areas.
- 5) Apply (with syringe without needle) a limited quantity of product on the skin in the area being treated.
- 6) Apply with the rubber tipped brush covered with a nitrile glove, repeating the application 2-3 times.
- 7) Proceed area by area until the entire desired surface has been treated.
- 8) Remove the residue with a pad dampened with water.
- 9) Restore the hydrolipidic film with WiQo Restoring Cream for External Delicate Areas.

In case of hyperpigmentation, continue the treatment with the home product WiQo Lightening Serum for External Delicate Areas.

REVITALISATION AND LIGHTENING OF MAMMARY AREOLA

Frequency of sessions:

1 session every 15 days for 5 or more sessions.

METHOD OF USE

- 1) Clean the area to be treated with WiQomed P Solution.
- 2) Perforate the cap of the vial with a needle (18g).
- 3) Insert the syringe on the needle and extract the desired quantity.
- 4) Apply on the areola with the rubber tipped brush covered with a nitrile glove, repeating the application 2-3 times (avoiding the nipple).
- 5) If irritation or intense burning occurs, stop the application.
- 6) Apply WiQo Lightening Serum for External Delicate Areas.

In case of hyperpigmentation, continue the treatment with the home product WiQo Lightening Serum for External Delicate Areas.

WiQO[®]med

REVITALISATION AND LIGHTENING OF AXILLA

Frequency of sessions:

1 session every 15 days for 5 or more sessions.

METHOD OF USE

- 1) Clean the area to be treated with WiQomed P Solution.
- 2) Perforate the cap of the vial with a needle (18g).
- 3) Insert the syringe on the needle and extract the desired quantity.
- 4) Apply with the rubber tipped brush covered with a nitrile glove, repeating the application 2-3 times.
- 5) If irritation or intense burning occurs, stop the application.
- 6) Apply WiQo Lightening Serum for External Delicate Areas and WiQo Restoring Cream for External Delicate Areas.

In case of hyperpigmentation, continue the treatment with the home product WiQo Lightening Serum for External Delicate Areas.

REVITALISATION AND LIGHTENING OF PERIANAL AREA

Frequency of sessions:

1 session every 7 days for 5 or more sessions.

METHOD OF USE

- 1) Clean the area to be treated with WiQomed P Solution.
- 2) Perforate the cap of the vial with a needle (18g).
- 3) Insert the syringe on the needle and extract the desired quantity.
- 4) Protect the mucous membranes with WiQo Restoring Cream for External Delicate Areas.
- 5) Apply with the rubber tipped brush covered with a nitrile glove, repeating the application 2-3 times.
- 6) If irritation or intense burning occurs, stop the application.
- 7) Apply WiQo Lightening Serum for External Delicate Areas and WiQo Restoring Cream for External Delicate Areas.

In case of hyperpigmentation, continue the treatment with the home product WiQo Lightening Serum for External Delicate Areas.

GENERAL CONTRAINDICATIONS IN THE APPLICATION OF PRX-T Lady

- Irritated skin
- Irritable skin
- Macerated skin
- Active dermatitis
- Allergy to kojic acid
- Allergy to latex, if latex gloves are used
- Pregnancy and breastfeeding
- Cutaneous and/or herpetic eruptions